NMR data processing software

**Delta**

NMR Software

v5.0

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
Experiment Parameters

Header Instrument **Acquisition** Pulse Diagram Favorites

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
The image displays two windows from a Bruker software interface. The main window, titled "Experiment Parameters", has a tabbed interface with "Header" selected. The "Header" tab contains several parameters: "storage\_filename" (test\_Proton), "filename" (proton), "storage\_comment" (single\_pulse), "comment" (single\_pulse), "auto\_gain" (checked), "filter\_limit" (16), "force\_dual\_mode" (unchecked), and "force\_tune" (checked). The "force\_tune" checkbox is circled in red. A smaller window titled "Job" is overlaid on the top right, showing a tree view of "test" (0h 02m) containing "Proton" (0h 02m). The "Proton" entry is circled in red with a "1" next to it. In the "Job" window's toolbar, an edit button (a document with a pencil) is circled in red with a "2" next to it. In the "Experiment Parameters" window, the "Header" tab is circled in red with a "3" next to it.

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The screenshot shows the 'Attributes' table for a sample named 'test'. The table has the following columns: Attribute Name, Value, and a checkbox. The 'gradient shim' attribute is unchecked, and the 'spin\_state' attribute is set to 'SPIN OFF'. Red arrows point to these settings with Chinese annotations: '不勾选即为不匀场' and 'Spin off为不旋转'.

Attribute Name	Value	Checkbox
gradient shim		<input type="checkbox"/>
last shimmed	4-DEC-2020 11:24:09	
load shims	System	
lock_state	AUTOLOCK	
preparation		<input checked="" type="checkbox"/>
spin_set	15[Hz]	
spin_state	SPIN OFF	<input checked="" type="checkbox"/>
temp_set	25[dC]	
temp_state	TEMP OFF	

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