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The screenshot displays the 1D Processor software interface. The top menu bar includes File, Options, Reports, PreTransform, Window, Transform, PostTransform, Display, Analyze, and Tools. The toolbar contains various icons, with the save icon (a floppy disk) circled in red. A red arrow points to this icon with the text "点击此处 保存处理后的数据" (Click here to save the processed data). The main window shows two plots: the top plot is a time-domain waveform with a decaying envelope, and the bottom plot is a frequency-domain spectrum with a peak at approximately 1.4 ppm. The bottom plot's x-axis is labeled "X : parts per Million : Proton" and the y-axis is "abundance". The file name "sample_Proton-1-2.jdf" is circled in red in the bottom plot's title bar. On the right side, there are two panels: "Processing Tools" and "Options". The "Processing Tools" panel shows a list of processing steps: ▶sexp(0.2[Hz], 0.0[s]), ▶trapezoid(0[%], 0[%], 80[%], 100[...], ▶zerofill(1, TRUE), ▶fft(1, TRUE, TRUE), machinephase, ppm, and [display/phase]. The "Options" panel shows parameters for Peak Sep. (0[Hz]), Int Width (50[Hz]), Normal (1), X Ref (0.0), X Start (12.50[ppm]), and X Stop (-2.50[ppm]).

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